

ImmunoStrip® for Bt-Cry1F, Bt-Cry2Ab2, Bt-Cry34Ab1 & Bt-Cry3Bb1



Part Number: STX 63200/0050

Test Label: NA

Test Format: ImmunoComb® or ImmunoStrip®

Capture Reagent: NA

Detection Reagent: NA

Intended Use:

This quad trait ImmunoStrip® test is intended for seed quality purposes to determine the presence or absence of insect resistance traits found in

SmartStax® corn seed. In just minutes, the test can independently detect Bt-Cry1F, Bt-Cry2Ab2, Bt-Cry34Ab1 and Bt-Cry3Bb1 proteins on one single strip when testing seed or leaf tissue.

This ImmunoStrip® test is offered in canisters of 50 strips or in 48-well and 96-well comb configurations for high throughput testing applications.

*SmartStax® multi-event technology developed by Dow AgroSciences and Monsanto.
SmartStax® and the SmartStax® Logo are trademarks of Monsanto Technology LLC.
ImmunoComb® and ImmunoStrip® is a registered trademark of Agdia, Inc.*

Includes:

- ImmunoComb® or ImmunoStrip®
- SEB7 sample buffer not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com