

ImmunoStrip® for 2mEPSPS



Part Number: STX 64100/0050

Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Monoclonal
Detection Reagent: Monoclonal

Intended Use:

This Agdia ImmunoStrip® for 2mEPSPS is designed and validated for detecting the presence or absence of the 2mEPSPS protein in single seed and leaves and for testing bulk mixtures. The Agdia 2mEPSPS ImmunoStrip® is validated for testing transgenic soybean with a detection limit of 1 in 200 seeds (0.5%).

These strips are commonly used by plant breeders, universities and private testing labs who want fast, clear results.

Please see product insert for more information.

Includes:

- ImmunoStrips®
- SEB4 sample buffer not included, see related products
- [User](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com