

ImmunoStrip® for CP4 EPSPS (RR)



Part Number: STX 74000/0050

This Agdia ImmunoStrip® for Roundup Ready® is designed and validated for detecting the presence or absence of the CP4 EPSPS (RR) protein in single seed and leaves and for testing bulk mixtures. The Agdia Roundup Ready® ImmunoStrip® is validated for testing transgenic corn, cotton, and soybean with a detection limit of 1 in 1000 seeds (0.1%).

These strips are commonly used by plant breeders, universities and private testing labs who want fast, clear results.

Roundup Ready® is a registered trademark of Bayer Group.

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com