

ImmunoStrip® for Vip3A



Part Number: STX 83500/0050

Test Label: NA

Test Format: ImmunoStrip®

Capture Reagent: Monoclonal

Detection Reagent: Monoclonal

Intended Use:

This ImmunoStrip® test is intended for seed quality purposes to determine the presence or absence of Vip3A protein in transgenic corn.

This test system can be used to test individual corn seed and corn leaf or detect 1 transgenic Vip3A seed in 500 corn seeds (0.2%).

The Vip3A ImmunoStrip® has shown no cross-reaction with other transgenic proteins in corn seed and leaf including CP4 EPSPS (Roundup Ready® corn), Bt-Cry1Ab, Bt-Cry1F, Bt-Cry2A, Bt-Cry34Ab1, Bt-Cry3Bb1, mBt-Cry3A, CspB, PAT/bar, GA21 and GT21.

Includes:

- ImmunoStrips®
- PBST sample buffer not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com