

## ELISA Reagent Set for Maize dwarf mosaic virus (MDMV)

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

---

This ELISA is a qualitative serological assay intended for *Maize dwarf mosaic virus*(MDMV) detection in plant foliage.

This test may have a mild cross-reaction with some isolates of *Sugarcane mosaic virus*. It does not cross react with *Johnsongrass mosaic virus*.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Maize dwarf mosaic virus (MDMV), 96 testwells,SRA 18000/0096

Reagent Set for Maize dwarf mosaic virus (MDMV), 500 testwells,SRA 18000/0500

Reagent Set for Maize dwarf mosaic virus (MDMV), 1000 testwells,SRA 18000/1000

Reagent Set for Maize dwarf mosaic virus (MDMV), 5000 testwells,SRA 18000/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)