

ELISA Reagent Set for Sugarcane mosaic virus (SCMV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:


This ELISA test is a qualitative serological assay intended for detection of *Sugarcane mosaic virus* (SCMV) in plant foliage.

The test can detect several SCMV strains, including USA A, B, D, E, F, Australian SCMV strains BC, Sabi, and SC.

This test may have a mild cross-reaction with some *Maize dwarf mosaic virus*. It does not cross-react with *Johnsongrass mosaic virus*.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Sugarcane mosaic virus (SCMV), 96 testwells, SRA 18100/0096

Reagent Set for Sugarcane mosaic virus (SCMV), 500 testwells, SRA 18100/0500

Reagent Set for Sugarcane mosaic virus (SCMV), 1000 testwells, SRA 18100/1000

Reagent Set for Sugarcane mosaic virus (SCMV), 5000 testwells, SRA 18100/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com