

## ELISA Reagent Set for Arabis mosaic virus (ArMV)

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

This ELISA is a qualitative serological assay for detection of *Arabis mosaic virus* (ArMV) in plant foliage.

The test is proven to detect over 33 ArMV isolates from all over the world.

The test has been evaluated by Agdia Testing Service Department in over 162 different plant hosts including ornamental, vegetable, and fruit species. Agdia ArMV ELISA test is used in our Testing Services Ornamental, Berry, Petunia, Hosta, and Rose crop screens.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- Tween 20
- General extract buffer 3, GEB3
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)

Reagent Set for Arabis mosaic virus (ArMV), 96 testwells, SRA 23203/0096

Reagent Set for Arabis mosaic virus (ArMV), 500 testwells, SRA 23203/0500

Reagent Set for Arabis mosaic virus (ArMV), 1000 testwells, SRA 23203/1000

Reagent Set for Arabis mosaic virus (ArMV), 5000 testwells, SRA 23203/5000

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)