

ELISA Reagent Set for BCMV

Test Label: Alkaline Phosphatase

Test Format: Antigen coated plate (ACP) ELISA

Capture Reagent: Polystyrene Plate

Detection Reagent: Monoclonal

Intended Use:

The intent of this test is to detect all strains of BCMV and BCMNV, generally the most important potyviruses affecting beans, worldwide. The test can be performed easily with any plant that can be infected with a potyvirus.

This is not a specific test for BCMV or BCMNV.

In order to detect all BCMV and BCMNV strains, the test also detects other potyviruses infecting beans.

Includes:

- Detection antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Bean common mosaic virus (BCMV), 96 testwells, SRA 46000/0096

Reagent Set for Bean common mosaic virus (BCMV), 500 testwells, SRA 46000/0500

Reagent Set for Bean common mosaic virus (BCMV), 1000 testwells, SRA 46000/1000

Reagent Set for Bean common mosaic virus (BCMV), 5000 testwells, SRA 46000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com