

## ELISA Reagent Set for Psph

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

**Synonym:** *Pseudomonas syringae* pv. *phaseolicola*

### Intended Use:

The ELISA is a reliable tool for *Pseudomonas savastanoi* pv. *phaseolicola*, the causal agent of halo blight of beans.

Symptoms include lesions on cotyledons, leaves, stems, and seed coats. Leaf lesions are usually surrounded by a yellow, chlorotic "halo" area. The test can be used with bacterial cultures, bean seed, or plant tissue.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 96 testwells, SRA 53300/0096

Reagent Set for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 500 testwells, SRA 53300/0500

Reagent Set for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 1000 testwells, SRA 53300/1000

Reagent Set for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 5000 testwells, SRA 53300/5000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)