

ELISA Reagent Set for WMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

Agdia WMV ELISA is a rapid, sensitive and economic serological assay for reliable detection of *Watermelon mosaic virus* (WMV) (formerly referred as WMV-2) in cucurbits (cucumber, melon, squash, pumpkins, and watermelon) and legumes (bean, pea, soybean, cowpea, etc).

The test does not cross react with other related viruses, including: AMV, CMV, CGMMV, CYYV, CYSDV, GRSV, KGMMV, MNSV, PRSV, SqMV, TCSV, ToRSV, TRSV, TSWV, WSMoV, and ZYMV.

Includes:

- Capture antibody
- Alkaline phosphatase conjugated antibodies
- Tween 20
- General extract buffer 3 (GEB3)
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Watermelon mosaic virus (WMV), 96 testwells, SRA 54001/0096

Reagent Set for Watermelon mosaic virus (WMV), 500 testwells, SRA 54001/0500

Reagent Set for Watermelon mosaic virus (WMV), 1000 testwells, SRA 54001/1000

Reagent Set for Watermelon mosaic virus (WMV), 5000 testwells, SRA 54001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com