



Transfer Guns

Also known as the IMV Mini Syringe.

F17 deep chamber accepts full length 0.25ml straws. F17A shallow chamber rod accepts 0.25ml straws that have been clipped or shortened immediately prior to insertion into the rod. Both rods require a 21" sheath, choose from the Agtech F20 or F21 or the IMV F18A or F18.

F17	Deep chamber
F17A	Shallow chamber
F23A	IMV Equine ET Syringe Gun.5ML (fits with F22 sheaths)

Agtech "Optimum" Sheath, 21"

- Stainless steel, side delivery tip sheath, preferred for transferring embryos.
- Fits the IMV transfer rod, items F17 and F17A.
- Fits F23A gun.



Comparable to the IMV sheath the Agtech embryo transfer sheath costs less, but performs equally as well. The sheath's stainless steel side-delivery port, through which the embryo is expelled, is less likely to become plugged with mucus and debris during the transfer procedure when compared to conventional transfer sheaths. The round smooth sheath tip minimizes uterine trauma during transfer.

F20	Pack of 5, no chemise. For Bovine or smaller embryos, 0.25ml straws
	1-29 packs
	30 or more packs
F21	Pack of 25, individually packaged w/chemise. For Bovine or smaller embryos, 0.25ml straws
	1-5 packs
	6 or more packs
F22	Pack of 5, individually packaged w/chemise. For Equine Embryos, 0.50ml side delivery sheath, 0.50ml straws
	1-5 packs
	6 or more packs



IMV Oversleeve Sanitary Chemise

Protective plastic sleeve is pulled over the transfer or insemination rod immediately prior to entering the cow/heifer. Prevents the transfer of contaminants into the uterus during insemination or embryo transfer. Upon contact with the cervix, the sleeve is pulled back forcing the clean rod through the sleeve.

F27B	18" oversleeve, roll of 100
F27C	21" oversleeve, roll of 80

Sani-Shield Rod Protectors

Sterile PVC "Sani-Shield Protector" is used as a protective cover for assembled rods. It prevents the transfer of contaminants into the uterus during insemination or embryo transfer. A slight forward push will enable the rod to pass through the protective tip. Individually wrapped, packed 25/poly bag.

B7-4525	8" (20.32cm) long
B7-4625	10" (25.40cm) long
B7-4725	12" (30.48cm) long

Insemination Gun, O-Ring

- Uses sheath with green insert
- For 0.25ml straws only

Made of durable stainless steel with welded components. Dimensions are 18" (45.72cm) long, .156" (.40cm) O.D. Use with sheaths B6-4400 or B6-4400C. Gun for 0.25ml straws.

7070	Gun for 0.25ml straws
------	-----------------------

Insemination Gun, Spiral Hub

- Uses spiral sheath with blue insert
- Universal gun accepts either 0.25ml or 0.5ml straws

Made of durable stainless steel with welded components. Dimensions are 18" (45.72cm) long. Use with sheaths B6-4500 or B6-4500D

7072	Universal gun
------	---------------

Insemination Gun & Sheath for Sheep & Goats

"All-2-Mate" universal gun accepts 0.50ml & 0.25ml straws.

G6-3250	AI gun
G6-4400	Sheath 25/pk

AI Sheaths with insert

- Insert provides a secure seat for 0.25ml and 0.5ml straws
- Dimensions 17.5" (44.45cm) long, .278" (.452cm)O.D.

Green insert sheath is slit and slides on the o-ring hub gun. Blue insert sheath is un-slit and easily twists on the spiral portion at the head of the spiral gun. 50 sheaths/pack.

B6-4400	Green insert sheaths, individually packaged, for O-Ring gun.
	4 or more packs
B6-4400C	Green insert sheaths, bulk package, for O-Ring gun.
	4 or more packs
B6-4500	Blue insert sheaths, bulk package, for O-Ring gun.
	4 or more packs
B6-4500D	Blue insert sheaths, individually packaged, for O-Ring gun.
	4 or more packs