



Vortech™ Silicone Bovine

Vortech™ Silicone Bovine Catheters 23" long, Sterile

- 100% pure medical grade silicone with 6 large off-set eyelets
- Sterile and autoclavable.
- Extra long length: 23" (58cm)
- Foley 2 way, closed end



Agtech Silicone catheters are sterilized for ready-to-use convenience. A stylet tip "pocket" at the end of the Agtech catheter surrounds and guards the stylet reducing chances of uterine trauma. Autoclavable and durable so you can use it again and again. A longer catheter simplifies horn flushing large, deep cows with heavy tracts. Available in 5cc or 30cc balloon. French sizes 14fr to 22fr.

E112	14fr 5cc balloon	E104	20fr 5cc balloon
E100	16fr 5cc balloon	E105	20fr 30cc balloon
E101	16fr 30cc balloon	E106	22fr 5cc balloon
E102	18fr 5cc balloon	E107	22fr 30cc balloon
E103	18fr 30cc balloon		

1-9 catheters of the same item number
10 or more of the same item number

Silicone Small Ruminant Catheter

Autoclavable, re-usable, and sterile silicone catheter for small ruminants.

ECE083	8fr 5cc 12" long
--------	------------------

Vortech™ Silicone Catheter Syringe Flush

- Sterile
- Foley 2 way, closed end
- Similar to the Minitube "Desi" catheter
- Used for Syringe Flushing



E203	18fr 30cc balloon, 25.2" (64.1cm)
	1-99 catheters
	100-499 catheters
	500+ catheters
E204	20fr 30cc balloon, 25.2" (64.1cm)
	1-99 catheters
	100-499 catheters
	500+ catheters
E205	18fr 40cc balloon, 33.1"(84cm)
	1-99
	100-499 catheters
	500+ catheters

Vortech™ Silicone Equine Catheters 34" long, Sterile

- 100% pure medical grade silicone with 4 large off-set eyelets, open end
- Sterile and autoclavable.
- Extra long length: 34" (86cm)



Silicone catheters from Agtech are sterilized for ready-to-use convenience. Autoclavable and durable so you can use it again and again. French sizes 28fr to 36fr.

E108	28fr 80cc balloon
E109	32fr 80cc balloon
E110	36fr 80cc balloon
	1-9 catheters of the same item number
	10 or more of the same item number



Latex Catheters 18" Long

Foley 2 way catheter is sterile & disposable. Silicone lined, latex catheter recommended one time use. Available in 5cc or 30cc balloon, French sizes 14fr (smallest diameter for use with E23 stylet) to 24fr (for use with stylets E24). An industry standard.

E02	14fr 5cc balloon	E05	20fr 5cc balloon
E09	14fr 30cc balloon	E12	20fr 30cc balloon
E03	16fr 5cc balloon	E06	22fr 5cc balloon
E10	16fr 30cc balloon	E13	22fr 30cc balloon
E04	18fr 5cc balloon	E07	24fr 5cc balloon
E11	18fr 30cc balloon	E14	24fr 30cc balloon
	1-9 catheters of the same item number		
	10 or more of the same item number		

Stylets

Stainless steel stylets for use with latex and silicone catheters.

E23A	25" long, small diameter, 14fr or smaller
E24A	25" long, large diameter, for 16fr or larger
E24AC	34.5" long, large diameter, for 16fr or larger

Stylette

Bioniche stainless steel safety coated stylette.

ECE020	23.6" long
--------	------------



E203



Vortech Silicone Equine

Latex 18"



E23A - E24A