

Filters are sterile, individually packaged and compatible with our Y-junction tubing, E27. Used for embryo collection during flushing.

Agtech zona™ Filter

- Tethered cap on the inflow port to securely seal post flush.
- Air vent with cap
- Pinch clamp for regulating outflow.
- 120ml capacity
- 75 micron stainless steel screen.
- Compare to the EmCon Filter.



D03	1-19 filters
	20-99
	case of 100

EmCon Filter

- 120ml capacity
- 75 micron stainless steel screen
- Includes drain tube with pinch clamp

D02	1-19 filters
	20 or more

Minitube/MiniFlush Filter

- Halves separate after collection to permit searching in the dish
- Includes drain tube

D05A	1-19 filters
	20 or more

VCI Filters

- High Volume 120ml media capacity
- Low Volume 50ml media capacity
- 75 micron mesh nylon filter

D04	High volume filter
	1-19 filters
	20 or more

D05	Low volume filter
	1-19 filters
	20 or more

Agtech Rapid-Flow™ Y-Tubing

- Radiated, individually packaged.
- Each 5' (152 cm) tubing arm includes a color coded pinch clamp.
- Spike end mates to 1 liter and 2 liter medium bags.
- Catheter end mates to bovine & equine catheters.
- Filter end mates with all filters listed.



E27	1-19 tubing
	20-99
	case of 100

Syringe Filter

Low Protein binding membrane maximizes recovery of critical components. Gamma radiated, individually blister packaged. Attaches to luer tip of syringe: 0.2 u filter screen, 25mm disc diameter. Manufactured by Starlab Scientific.

D06G	Starlab Scientific	10 filter pack
D06F	Starlab Scientific	100 filter pack

IVF (In Vitro Fertilization) Catheter

- For aspirating and manipulating embryos between dishes.
- Packaged 5 radiated catheters per individual zip lock bag.
Sold 50 catheters /pack. Use w/ item G11, 1ml syringe.

D17	1-5 packs
	6 or more packs

Tom Cat Catheters with Suture Adapter

Individually packaged & sterile. 50 catheters/box.

703013	4.5" long open end. 1-9 boxes
	10 or more boxes
703021	5.5" long open end. 1-9 boxes
	10 or more boxes

Mechanical Pipettor

Lightweight, fixed volume 5ul full displacement pipettor now with a built in tip ejector.

EHP11	5ul Pipettor
EHP12	1-200ul capacity tips. Rack of 96 sterile tips
EHP12A	Snowflake Style

Microdispenser, fixed volume

Delivers fixed volume at increments of 1,2,3,4, or 5ul using glass capillary bores.

D22	Model D22 microdispenser
D23	Glass bores for model D22. 100/pk, sterile
D21	Glass bores for model D20. 100/pk, sterile

