



**B4-7430 - B4-7609**



**F30 - F31**



**F09A - F09AY**



**TLS-2200**



**F10**



**C19391N**

**61625**



**SAVB-104**



**SAVL-102-25**



**C15245N**



**SAVF-102-25**

**Goblets**

Durable plastic goblets attach to canes for holding straws in nitrogen tanks. 10mm goblets hold 3-4 .5cc straws or 4-5 .25cc straws. 13mm goblets hold 5-6 .5cc straws or 6-7 .25cc straws.

B4-7430	10mm clear goblets. 200/pack
B4-7439	10mm yellow goblets. 200/pack
B4-7440	13mm clear goblets. 150/pack
B4-7609	13mm yellow goblets. 150/pack

**Canes, Aluminum**

Aluminum canes for holding goblets.

F30	11" cane holds two 10mm goblets, 100/pk.
F30C	11" cane holds two 10mm goblets, Case of 1,000
F31	11" cane holds two 13mm goblets, 75/pk.
F31C	11" cane holds two 13mm goblets, Case of 1,000

**Cane ID Tabs**

Aluminum tabs for marking 10mm or 13mm canes. Pack of 100 tabs.

F09A	white
F09AY	yellow

**Thermal Labeling System**

An easy and legible method for labeling embryo straws, goblets and canes. The hand-held, battery operated TLS-2200 allows you to quickly and effortlessly create sharp professional labels for:

- .25cc straws
- .50cc straws
- Sealing plugs
- Plastic goblets and aluminum canes

Label material, adhesive and ink are 100% stable in liquid nitrogen and will not detach or rub off. Clearly a superior method for labeling your valuable embryo containers. Weighing in at just 2.75 lbs (1.25kg) it's the only bar code/label printer that you can hold in just one hand.

Unit is powered by a detachable 12v battery which is rechargeable when connected to a 110v electrical source.

TLS-2200 labeling system includes: Labeler, printer ribbon, rechargeable 12v battery (use 110v for charging.)

TLS-2200

**Labels and parts for TLS-2200 Labeler**

Thermal labels are 100% nitrogen stable. Each roll contains 250 labels.

PTL-19-427	White labels for .25cc straws. 1" x .5" print area
PTL-19-427-YL	Yellow labels for .25cc straws. 1" x .5" print area
PTL-32-427	White labels for .25cc straws. 1.5" x .5" print area
PTL-32-427-YL	Yellow labels for .25cc straws. 1.5" x .5" print area
PTL-1375	Yellow labels for .25cc straws, 1.375" x 1.0" print area
PTL-19-642	White labels for 0.5cc semen straws, 1" x .5" print area
R4310	Replacement ribbon, black
PCK-5	Printer head cleaner kit

**Marking Pen**

Ultra fine, permanent black ink, "Sharpie" marker.

F10

**AV Cover, Neoprene**

Insulated cover with four Velcro straps. Washable and holds AV better than leather.

C19391N



**16" Missouri A.V.**

- Lighter and easier to handle than the CSU model

61625



**16" Missouri Model A.V. Kit**

- New filling hose, cotton filters

Comes complete with 16" latex AV, neoprene cover, two collection bottles, bottle adapter, cotton semen filters, no spray filling system.

61629



**Colorado Bladder**

- Molded from natural latex rubber
- Acts as a water retention barrier inside the AV Case

Two 7" x 5/8" rubber bands are required to secure the bladder to the AV case.

SAVB-104 22" flared end bladder



**Colorado Disposable Liners**

22" Soft, flexible plastic film with tapered end to connect to the collection bottle and a flared end for easy assembly. Extra soft plastic provides a softer feel, while standard provides superior durability. From ARS.

SAVL-102-25	Extra soft, 25 liners
SAVL-101-25	Standard, 25 liners



**Semen Cotton Filter**

Equine artificial vagina semen cotton filter - 6" x 2-1/4"

C15245N 25/bag

**Semen Filters, Equine**

- Non-absorbent nylon mesh
- Fits inside disposable collection bottle

This ARS filter allows the sperm rich fraction of the ejaculate to pass into the collection bottle while the excess gel, smegma and debris are captured for easy disposal.

SAVF-102-25 pack of 25