



586A-201



B6-6280



C52



H05E



E176



H31



J18

J20

***Water Bath**

Excellent for thawing frozen semen straws, and warming the E-Z Mixin semen extender, or any other equipment that requires a fixed physiological temperature. The thermoregulated water bath, from ARS, is factory calibrated to maintain water temperature at 37C +/-0.5C (98.6 +/-0.9F). Its convenient size (12.8"L x 11.5"W x 7"H) takes up minimal counter space, making it ideal for laboratories and indoor work areas. The water bath kit includes an instructional manual, transparent plastic lid and a detachable power cord.

586A-201 Water Bath Kit

Thermos with Thermometer

An unbreakable, light weight, rugged, wide mouth plastic thermos bottle with molded side handle and a capacity of one pint (671.30ml). Includes a thermometer. Dimensions are 4.5" (11.43cm) wide and 8.0" (20.32cm) long.

B6-6270 Celsius

B6-6280 Fahrenheit

Thaw Unit, IFT

A simple to use and user friendly semen thaw unit from IFT, which is built to survive under any working conditions. Robust unit is made out of high impact plastic and has no glass inner flask. Automatically heats and maintains an accurate temperature for semen thawing. Pre-Set temperature at 36°C.

The Inner cylinder contains a 4 slot divider which helps you identify/organize your semen straws. The Lifter elevates the straws in order to prevent the use of fingers which can cause the semen to warm unevenly, lowering the number of live sperm.

- **No need to worry about voltage source:** 3 in 1 - 12V/110/220V (power pack kit) 4 meter power cord / 12V car plug adapter. 12V socket with battery connecting clips. 110/220V automatic power supply wall adaptor for mains.
- **Easy to See - Solid Stem Thermometer** for rapid reading of temperatures. Shockproof protection sleeve cover, -40°C / 50°C - 16cm.
- **No Inner Glass flask**
- **Drop Proof**
- **4 Slot chamber**
- **Straw Lifter**
- **CTT™ heating system**
- **Carry case**
- **2 year warranty**

C52 Thaw Unit

Thermocouple Thermometer

Temperature can be displayed in degrees F or C as large, 4-digit LCD numbers. HOLD feature lets you freeze the display for later reading. Thermometer is calibrated and conforms to NIST tables for accuracy. Probe is ideal for verifying temperatures in your embryo freezer or warm-water bath. Range: -190 to 1000c; Resolution: 0.1 c; Accuracy: +/- 0.1 c; Display: 4-digit LCD, Power: 9V battery/1000 hours; Probe: stainless steel, 4.5" Long, 1/8" diameter, with cable; Dimemnsions: 3.3" W x 6.3" H x 1.1" D; Warranty: 3 years

H05E

Thermometer, Digital

Digital Thermometer, Temp. Range 14F (10C) to 392F (200C). Dial Size 1 3/4 In., Dia. 1 3/4 In., Length 6 In., Stem Length 5 In. Construction Plastic, Includes Battery. Recommended for use with the C52 Thaw Unit.

B6-6330D

Thermometer

E176 Extra Long-Stem Thermometer 14 3/8", digital C & F, stainless steel probe

ET & AI Supply Case

What a great way to keep track of your embryo transfer and a.i. supplies. This molded heavy-duty "structolene polyethylene" case is impervious to chemicals, manure and moisture. Double heavy-duty plated latches reliable holds case closed. The case contains a compartment tray to hold small items like your straws scissors and measuring tape. 26" long x 8 3/4" wide x 10.5" high.

H31

Sterilization tubing

This sterilization packaging material is clear plastic on one side and paper on the other. Cut to length and either heat-seal or tape closed. 100' roll.

J18 3' wide roll

J20 6' wide roll



Estroprotect™ Heat Detector Patch

Excellent visual aid for detecting cows in standing heat. Applies by self-adhesive to the tailhead. When cows are mounted, the scratch off surface is rubbed off exposing the brilliant signal layer below. 50 markers/pack.

J13	Orange
	2 or more packs
J14	Green
	2 or more packs

Kamar Heat Detector Patch

Kamar heatmount detectors allow you to identify more cows in heat and to create more profitable calving intervals. To use the Kamar heatmount detector, glue the detector onto the tailhead of a cow. When a cow wearing a Kamar is mounted by a herd mate constant pressure from the brisket of the mounting animal turns the detector red. This leaves a visible indication that the cow stood to be ridden and therefore may be in heat and ready to breed. An exclusive built-in timing mechanism is designed to distinguish false mounts from true standing mounts. Box of 50 detectors includes 2 tubes of adhesive.

J15	1 box
	2 or more boxes

Sharp Mark marking crayons

All weather marking crayon for hipmarking recipients, culls and treated cattle. 12 crayons/box.

J24C	Fluorescent Orange
------	--------------------

Embryo Data Form

Organize all your important embryo collection, transfer and freezing information, each time you flush on this single form. Format approved by the International Embryo Transfer Society.

V26	20 forms/pack
-----	---------------

Kimwipes

Kimwipes are an extra low linting tissue that will not scratch or contaminate most delicate surfaces. They will easily wipe up both liquid and dust and are highly absorbent. The dispenser is anti-static to reduce electrostatic charge. Each box contains 280 wipes.

06-666A	1 box
	10 or more boxes

Speculum, Equine

The disposable vaginal speculum is 18" long with foil lining for greater visibility. The penlight illuminates the interior of the speculum for improved visibility while performing A.I. or other exams. Disposable vaginal speculums are sterile, individually wrapped and sold in packages of 10.



EDVS-102-10	10 speculums
537-664	speculum penlight

Culture Swab

A unique tool for collecting uterine samples from brood mares for cultures. Three piece design includes: an outer clear PVC tube, a blue PVC swab sheath inside the outer tube, a white PVC swab holder with attached cotton swab. Sterile, individually wrapped.



E9-5200	25 swabs/bag
---------	--------------

Dispensing Bag

This dispensing bag has an easy open and easy close flip cap fill port. Hangs easily in your barn or clinic. Dose not include tubing. 1500 ml capacity. Perfect for use with C06, D-PBS Modified or C09A Powdered D-PBS, when making a flushing solution.

E28	1 bag
-----	-------

Forceps

Stainless steel forceps or hemostats. 5.5" long.

H23	Straight
H24	Curved

Scissors

5.5" stainless steel scissor has a blunt/sharp end.

H22	
-----	--

Alconox, detergent

For cleaning lab ware and instruments. Biodegradable, rinses clean.

J23	4lb carton
-----	------------

Nolvasan

Nolvasan (chlorhexidine diacetate), the only EPA-registered chlorhexidine disinfectant, works against at least 60 different bacteria, fungi, yeasts, and viruses. Nolvasan is non-corrosive, has minimal to no skin irritation, and retains antimicrobial activity in the presence of organic matter. Nolvasan's unique binding to skin proteins provides residual activity for as long as two days. Usual dilution is 1 oz per gallon of water.

J26	1 gallon
	4 or more

Chlorhex

A virucidal, bactericidal, fungicidal sanitizing solution useful in the disinfection of equipment and premises. Usual dilution is 1 oz per gallon of water. Chlorhexidine Gluconate.

J26A	1 gallon
	4 or more, gallons



J13

