



* Catalog Items with a * are not stocked. Please allow 7-9 days for product delivery.



HO8A



***Meiji EMZ-5**

- Recommended for embryo work.
- Continuous zoom magnification from 15x to 90x.
- 20x Eyepieces
- Image remains clear & sharp throughout entire range
- Crisp, high resolution, precision optics provide an excellent, erect, un-reversed stereoscopic image.
- Rugged, all metal body designed to last a lifetime.
- Binocular eyetubes inclined at 45 degree for fatigue free viewing.
- Reflected LED light cool stage.
- Lifetime Warranty

HO8A	EMZ-5
HO8AB	spare bulb, 5v 30w for non LED illuminators
H32	Additional Hard Sided Carrying Case

***Leica S6E Stereo Zoom**

- Recommended for embryo work.
- 6.3:1 zoom and 38° viewing angle
- 20x eyepieces provide 36.5mm field of view
- Cool, detachable fiber optic illumination
- 5 year Warranty

S6E	S6E microscope
H08BB	spare bulb



S6E

***Meiji EMZ-8TR Trinocular Stereo Zoom**

Trinocular stereo zoom microscope for embryo evaluation and video or digital image documentation. Allows for group viewing when using the optional MEIJI VDEMZ video kit and when connected to any standard monitor. Digital image capture is simple. Continuous zoom magnification from 15-90x w/20x eyepieces. Reflected LED light. Lifetime warranty.

HO8TR	Trinocular Microscope
-------	-----------------------

***Meiji Phase Contrast**

Compound microscope with phase contrast for semen evaluation and video or digital image documentation. Features a large flat top mechanical stage, objectives of 10,20, 40 & 100x, 6v 30W illuminator, 10x eyepieces. Lifetime warranty.

ML2870	Binocular Phase Contrast
ML2970	Trinocular Phase Contrast



HO8TR

***MICRA Microscope, Handycop**

- High-Quality Professional Lenses
- Portable Design and Simple Operation
- Bright and Long-Life LED Lighting
- Long-Life Battery
- Specially Designed Polycarbonate Slide

The A.I. Handycop is a high-powered portable microscope designed as a fast and reliable test kit for fertility screening in artificial insemination applications. Using the A.I. Handycop, the most active and healthiest sperm samples can be selected for efficient artificial insemination resulting in a high rate of successful pregnancy. The instrument is easy to handle and use, and the high-powered LED light guarantees clear and bright images even in dark working environments. It is powered by a lithium battery and can be carried anywhere. 400x.

A4000	
-------	--



A4000

***Meiji Video & Digital Photography Kit**

Transmit live color images to a monitor (not included) for group viewing and instruction, or record on tape. Or attach your Nikon Coolpix, Olympus Camedia, Canon Powershot, or Fuji Finepix digital camera to instantly record images for emailing or archiving. Kit includes color video camera, cables, video C-mount adapter and camera adapter. Please specify your digital camera model when ordering.

VDEMZ	For the EMZ-8TR stereo zoom microscope
VDML	For the ML5955 compound microscope.


Glass Slides and Cover Slips

Glass slides are 25mm x 75mm x 1mm.

L15C	80 slides
L16C	100 cover slips



VDEMZ



L15C

L16C

