

Infusette™ Tubes, Bovine

Infusette™ tubes/pipettes with drilled syringe end for easy coupling with a luer tip syringe and eliminates the need for an adapter. 25 tubes/pk.

B8-5230	21" non-sterile
---------	-----------------

Infusion Pipette, equine blue flexi-syringe connection

These pipettes have an attached syringe adapter and are specially designed for use with many of the common uterine infusion products. 25 pipettes/pk

E6-5315	22" sterile indiv. wrapped
E6-5415	25" sterile, indiv. wrapped


Semen Infusion Pipette, flexible tip, equine

- Rounded and polished entry tip to minimize tissue trauma
- Individually package and radiated.
- 25 pipettes/pack

E057	25"
------	-----


Flexible Insemination Tube

Appropriate for IUI techniques. 10 Tubes.

537-682-10	30"
------------	-----


Universal Pipette, Equine

- Suitable for fresh and frozen semen.
- Smooth round flex tip for gentle insertion.
- Precision tooled tip with inner cone for straws.
- Flexible to adapt to cervical conformation.
- For intra-uterine insemination (IUI).
- Frozen Semen: Use with 0.5ml straws and stylet. Pipette remains in place in the mare when multiple straws are required for insemination. Remove empty straw with the stylet and reload in place.
- Fresh/cooled Semen: Fill syringe with semen and connect syringe to Universal Pipette

E1275	75cm long, not sterile w/internal pipette
-------	---


Gloves, Exam

These latex wrist length gloves are powdered so they go on easy. Not sterile. 100 gloves per box.

J07	Medium
	2 or more boxes
J08	Large
	2 or more boxes

Sleeves, Shoulder Length

A 34" (86.36cm) long, 1 mil thick, disposable shoulder length glove manufactured of high quality polyethylene. A sturdy, airtight, universal size glove packed in an easy dispensing hard box.

B6-1010	100 sleeves/box
	2 or more boxes

Sleeves, Super Sensitive™

A 34.75" long shoulder length glove manufactured of high quality 1 mil soft light blue polyethylene. Designed to fit the hand with minimum material overlap, and has the ability to stretch providing exceptional tactile sensibility.

V7-2312	Packed 100/ dispenser box.
	10 or more boxes
J02A	10/bag, radiated, individually wrapped, sterile
	10 or more bags

Sleeves, Fingerless

Fingerless sleeve protects your arm while allowing you to feel with your fingers. Shoulder length sleeves are 1.25 mil thick. 100 sleeves/bag.

J04	10 or more bags
-----	-----------------

Lubogel-V

A lubricating jelly composed of a solution of Sodium Carboxy MethylCellulose in Propyleneglycol. Minimizes irritation to the skin and the delicate mucous membranes of animals. Not sterile.

6655	1 gallon
	4 or more, gallons

Therio-gel® Veterinary Fertility Lubricant

Lubricant for all animal species clinically shown to not harm sperm, embryos or sensitive tissue.

Therio-gel is perfect for:

- Artificial insemination of cattle
- Collecting semen from bulls
- Artificial insemination of mares
- Collecting semen from stallions
- Breeding soundness evaluations of stallions
- Inseminating sows, sheep, goats, llamas and alpacas
- Collecting semen for breeding pigs, sheep, goats, llamas and alpacas

6142	3.6oz tube (106 ml; 102 g)
	12 or more tubes
6142A	4 single-use packets (0.24oz / each packet)


B8-5230

E6-5315
E6-5415

E057

E1275
J07-J08

B6-1010

V7-2312
J02A
J04
