

**Straws with Wick & Powder (W&P), Radiated**

133mm gamma radiated IMV straws with no ethylene oxide risk. The RWB (red/white/blue) W&P straws have a plug that seals stronger when wetted properly, lowering the potential for "blow outs." Straws are packaged 5/ziplock bag; 10 bags (50 straws) /pack.

F01	.25 cc clear straws 3 or more packs
F01Y	.25cc yellow straws,for freezing DT embryos 3 or more packs
F03	.25cc clear RWB straws,for freezing DT embryos. Use w/yellow plugs. 3 or more packs
F02	.5cc clear straws 3 or more packs

**Straws with Wick & Powder (W&P), Radiated, DT**

Yellow .25cc IMV straws, sterile, printed DT for direct transfer. 5 straws per bag; 4 bags (20 straws) per pack. Special order item.

DT6300	.25cc yellow, DT 10 or more
--------	--------------------------------

**Straws W&P, Bulk, Not Radiated**

Standard wick and powder seals one end of this IMV straw. Package of 2000 straws, not radiated.

F04	.25cc clear RWB straws.
F06	.25cc clear straws
F04Y	.25cc yellow straws
F07	.5cc clear straws

**Agtech straw Connector**

Yellow straw connector from Agtech is versatile allowing for one end to fit inside a .25cc straw and the other end to fit inside a .5cc straw. Connectors are 18mm long.

Connectors are packaged 5 connectors per zip-lock bag and then are radiated; 10 bags/50 connectors.

F15B	Agtech adaptable connector 18mm long, yellow 3 or more packs
------	---

**Agtech straw Sealing Plugs**

Packaged 5 connectors per zip lock bag, radiated; 10 bags (50 connectors) per pack.

F15WT	Agtech .25cc plug. 47mm long. White 3 or more packs
F15LL	Agtech .25cc plug. 47mm long. Yellow 3 or more packs



**Sealing Powder**

Instantly seals the non-W&P end of straws when it comes into contact with water. Powder is polyvinyl chloride. 1 lb. bag.

B4-1050	
---------	--

**Tweezers**

A correctly designed 6.5" (16.51cm) long polyester tweezer for quick and easy removal of straws and goblets from canisters. They are shaped with a specially grooved tip that grips and surrounds the straw, yet does not crush it. Plastic transmits little heat or cold. Individually wrapped.

F14A	for .5cc straws
F14B	for .25cc straws

**Straw Cutter**

The Cito straw cutter is designed to quickly cut off the sealed end of straws. After insertion of straw, a push of a button cuts the straw at the proper length. Disassembles for cleaning and disinfection. It is made of non-corrosive plastic with a stainless steel blade and is 2.25" (5.72cm) in diameter.

F16	
-----	--

**Heat Sealer, Electric**

8" sealer used for thermal sealing of straw ends and sterilized tubing.

F11	110 volt
F11A	220 volt

**Repair Parts**

F12	Top Teflon cloth
F12A	Metal heating element
F12B	Insulating underpad

**Equine Vitrification Kit, Bioniche**

- First commercial equine vitrification kit in the world!
- Freeze equine embryos without the need for additional electronic equipment.
- Easy to use kit comes complete with detailed instructions



Formulation : One x 8ml vial of each: VS1-Vitrification Sol'n,VS2-Vitrification Sol'n, VS3-Vitrification Sol'n, VS Diluent-Vitrification Diluent. Note: The manufacturer highly recommends the use of RWB 1/4cc straws, item F03 page 15, because of their resistance to splitting during the freezing process. No refrigeration required.

EQU030	
--------	--



F04 - F07



F15B



F15WT



F15LL



B4-1050



F14A - F14B



F16



F11 - F11A



F12 - F12B



EQU030